

# RIDA®QUICK Norovirus (N 1403)


**Fig. 1**

**Fig. 2**

**Fig. 3**

**Literature**

1. Mead PS, et al. Food- related illness and death in the United States. *Emerg Infect Dis* 1999; 5:607-625.
2. Glass RJ, et al. The epidemiology of enteric caliciviruses from humans: a reassessment using new diagnostics. *J Infect Dis* 2000; 181(Suppl 2):S254-261.
3. Kaplan JE, et al. An outbreak of acute nonbacterial gastroenteritis in a nursing home: demonstration of person-to-person transmission by temporal clustering of cases. *Am J Epidemiol* 1982; 116:940-948.
4. Johnston CP, et al. Outbreak management and implications of a nosocomial Norovirus outbreak. *CID* 2007; 45:534-540
5. Corwin AL, et al. Shipboard impact of a probable Norwalk virus outbreak from coastal Japan. *Am J Trop Med Hyg* 1999; 61(6):898-903.
6. Evan HS. General outbreaks of infectious intestinal disease in England and Wales: 1995 and 1996. *Commun Dis Public Health* 1:165-171.
7. Chan MCW, et al. Fecal viral load and Norovirus associated gastroenteritis. *EID* 2006; 12: No. 8

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# RIDA®QUICK Norovirus (N 1403)

## 1. Intended Use

For *in vitro* diagnostic use. The RIDA®QUICK Norovirus Test is a quick qualitative immunochromatographic test for determining Genogroup 1 (GG I) and Genogroup 2 (GG II) noroviruses in stool samples. It is used as an aid in the diagnosis of gastroenteritis and for the analysis of

stool samples from children and adults with the symptoms of suspected gastroenteritis caused by the Norovirus.

## 2. Explanation

Noroviruses are a major cause of gastroenteritis worldwide with an estimated 23 million cases a year in the USA (1, 2). They are frequently involved in outbreaks in communal facilities, such as nursing homes, hospitals, day nurseries and prisons, and on cruise ships (3, 4, 5). Outbreaks caused by Noroviruses are reported more often than outbreaks caused by bacterial pathogens and they have a considerable impact on public health (6). The RIDA®QUICK Norovirus Test, which is based on

monoclonal antibodies, makes it possible to determine antigens in stool samples quickly and reliably and therefore speeds up patient management. The quick test is a straightforward, sensitive method for determining the antigens of both Genotypes I and II of the Norovirus. It is particularly suitable for use on small sample series.

## 3. Test principle

This quick test is a multi-step immunochromatographic membrane test which uses specific monoclonal antibodies directed against Noroviruses. Beside the sample well, the test cassette has a reaction window within which two horizontal lines are placed with immobilised antibodies. The test line (T) contains antibodies against Norovirus antigens. The control line C contains anti-mouse IgG antibodies. Conjugate 1 contains biotinylated antibodies against Norovirus antigens and Conjugate 2 consists of streptavidin with horseradish peroxidase bonded to it. For the test procedure, a defined quantity of the supernatant from a pre-prepared sample suspension is mixed with a defined quantity of Conjugate 1 and pipetted into the second smaller window of the test cassette (sample well). During a 10 minute incubation time at room temperature, the sample/conjugate mixture migrates from the filter pad which has absorbed the sample via the test membrane. While doing so, the antigen-conjugate

complexes in Norovirus-positive samples bond to the immobilised anti-Norovirus antibodies on the test line and the biotinylated antibodies with no antigens bond to the control line. After Conjugate 2 has been added and an incubation time of 1 minute at room temperature, the streptavidin attached to the peroxidase bonds to the biotin which has been immobilised via the specific antibodies.

The non-bonded peroxidase is then removed by adding Wash buffer to the reaction window. If the test is Norovirus positive, a blue line appears on the test line (T) within three minutes of adding substrate, along with a blue line on the control line (C). If the control line does not turn blue, this indicates that test has not worked properly and the result cannot be evaluated.

## 4. Materials provided

There are enough reagents in the pack for 20 determinations

Cassette	20 det.	20 single pouched test cassettes
Diluent	30 ml	sample dilution buffer, protein-buffered NaCl solution; contains 0.05 % Kathon CG; ready for use; blue coloured
Wash	10 ml	Wash buffer, phosphate-buffered NaCl solution; ready for use
Conjugate 1	7 ml	Biotin conjugate antibodies against Norovirus in stabilised protein solution; ready for use; coloured blue
Conjugate 2	5 ml	Streptavidin peroxidase conjugate in stabilised protein solution; ready for use
Substrate	7 ml	Hydrogen peroxide/TMB; ready for use
Pipet	25 pieces	Bag containing 25 disposable graduated pipettes

## Explanation of symbols

 IVD	In vitro diagnostic reagent	 REF	Article number
	Follow the instructions		Number of tests
 LOT	Lot number		Date of manufacture
	Expiration date		Manufacturer
	Store at		

## 5. Materials required but not provided

- small test tubes for stool suspensions (e.g. Eppendorf reaction tubes) with rack
- Vortex mixer
- Disposable gloves
- Stop clock
- Waste container containing 0.5 % sodium hypochlorite solution
- Disposable inoculation loops or disposable spatulas for stool sampling

## 6. Precautions

1. For *in vitro* diagnostic use only.
2. This test must only be carried out by trained laboratory personnel. The guidelines for working in medical laboratories must be followed.
3. The instructions for carrying out the test must be strictly adhered to.
4. The sample dilution buffer contains Kathon CG as a preservative. This substance must not be allowed to come into contact with the skin or mucous membrane.
5. Do not pipette samples or reagents by mouth.
6. Avoid contact with injured skin or mucous membranes.
7. When handling the specimens, wear disposable gloves and when the test is finished, wash your hands.
8. Do not smoke, eat or drink in areas where samples are being used.
9. All reagents and materials which come into contact with potentially infectious samples must be treated with suitable disinfectants (e.g. sodium hypochlorite) in exactly the same way as the samples themselves or autoclaved for at least one hour at 121°C.
10. The reagents in the RIDA®QUICK Norovirus Kit must not be used beyond the use-by date. All the reagents are adjusted for optimum activity. Diluting or changing these reagents will reduce their activity. Not keeping to the specified times and temperatures when carrying out the test will lead to false results.
11. Reagents from different lots of the RIDA®QUICK Norovirus Kit must not be mixed and used.
12. The reagents and samples must be protected from microbial contamination otherwise the results may be compromised.
13. A separate pipette must be used for each sample to prevent cross contamination and false results.
14. Once used, the cassettes must not be used again.

## 7. Stability and storage

The reagents must be stored at 2 – 8°C and must be used before the printed expiry date. After the expiry date, the quality guarantee is no longer valid. Likewise, we cannot guarantee that the cassettes are fit for use if

the foil packaging of the individual cassettes has been damaged so that moisture can penetrate.

# RIDA®QUICK Norovirus (N 1403)

## 8. Indications of instability or reagent expiry

The test must only be evaluated if the test cassette is intact **before** the sample suspension is pipetted into it and, except for a weak red colouration on the C side, no colour changes or bands are visible on the membrane. Furthermore, at least the **blue** control band (C) must be visible **after** the test incubation. If this does not appear, the following points must be checked before repeating the test:

- Has the expiry date on the test cassettes or the sample dilution buffer out of date?
- Has the correct test procedure been used?
- Is the extraction buffer contaminated?

If the control band is still not visible after repeating the test with a new test cassette, please contact the manufacturer or your local R-Biopharm distributor.

## 9. Sample collection and handling

Stool samples must be collected in clean containers without any additives and stored at 2 – 8°C before beginning the test. Sampling should take place as soon as possible after the symptoms appear (diarrhoea and vomiting) since virus elimination is at its maximum 1 – 3 days after the onset of the symptoms. If the sample is stored for more than 3 days, it must be frozen at – 20°C. Before starting the test, frozen samples must be

thoroughly thawed, brought to room temperature and (like the fresh samples) must be well mixed before use. Avoid freezing and thawing the specimen repeatedly. If rectal swabs have to be used, make sure that sufficient stool material (approx. 100 mg) is collected to carry out the test.

## 10. Preparation of samples and reagents

Before use, the samples, sample dilution buffer, reagents and test cassettes must be brought to room temperature (20 – 25°C). The test cassettes must only be removed from the external packaging shortly before

they are used. Once the cassettes have been removed from the packaging or used, they must not be used again. The test must not be carried out in direct sunlight.

### 10.1. Preparation of sample suspension

Place 1 ml sample dilution buffer [Diluent] in a labelled test tube. With the **liquid** stool sample, pipette 100 µl (up to just above the second thickening) of the sample with the disposable pipette [Pipet] and suspend it in the sample dilution buffer placed in the tube beforehand. With **solid** stool samples, suspend 50-100 mg (size of a small pea) in the buffer. The sample must then be well

homogenised. This can be achieved either by repeated suction and ejection of the suspension using the disposable pipette [Pipet] or, alternatively, by mixing it on a vortex mixer. After this, leave the homogenised suspension to settle for at least **2 minutes**.

## 11. Test procedure

**Figure 1:** picture of the test cassette and the disposable pipette

1. Remove the test cassette [Cassette] from the external packing and lay it on a level mat.
2. Remove 250 µl of the stool suspension supernatant using the disposable pipette already used for this sample (first mark) and place it in a clean, labelled test tube.
3. Add 6 drops (holding the vial vertically) [Conjugate 1]
4. Thoroughly mix the sample and pipette it complete with a slow steady stream into the sample well of the test cassette using the disposable pipette [Pipet]. Guide the pipette at an angle of approx. 45° with the point towards the reaction window so that the mixture can be placed fully on the membrane of the reaction window (Fig. 2).

**Figure 2:** Pipetting in the sample using the disposable pipette

5. Incubate the cassette for **10 minutes** at room temperature (20 – 25°C). Make sure that the membrane has been thoroughly saturated with the sample mixture. Otherwise, you will need to pipette a further 100 µl **[Diluent]** into the sample well.
6. Add 4 drops (holding the vial vertically) **[Conjugate 2]** to the reaction window and incubate the cassette for 1 more minute at room temperature.
7. Add 10 drops of wash buffer **[Wash]** to the reaction window and wait until the wash buffer is fully absorbed.
8. Add 6 drops (holding the vial vertical) substrate **[Substrate]** to the reaction window.
9. Read off the result within the next 3 minutes. Colourations of bands which appear after three minutes have no diagnostic value.

## 12. Interpretation of results

1. The interpretation of the test results is most reliable when the results are read off within the first 3 minutes of adding the substrate in good light and with a direct line of sight.
2. Check whether a blue line (the so-called internal control line) is visible on the “C” side of the reaction window and whether a blue line (the so-called test line) is also visible on the “T” side of the reaction window. The colour intensity of the lines may vary from weak to strong.
3. **Positive result:** A positive result can be recognised within 3 minutes of adding the substrate by the fact that both blue lines are visible (the test line “T” and the control line “C”). The colour intensity can vary from weak to strong. A clear part-line must be interpreted as positive. Other membrane colourations or shades on the other hand must **not** be interpreted as positive
4. **Negative result:** After 3 minutes have elapsed, a test can be interpreted as negative or not valid.

A single blue line is visible on the control side “C” of the reaction window but no blue line is visible on the test side “T” of the reaction window. A negative result indicates that there are no Norovirus antigens in the sample or the quantity of antigens is so small that it falls below the detection limit of the test.

5. **Result not valid:** There is no line on the test side “T” nor on the control side “C” of the reaction window or there is only a single line on the “T” side. In both cases, the test is invalid and should be repeated with a new cassette. The different interpretation models are illustrated in Figure 3.

**Figure 3:** Interpretation of results from the RIDA®QUICK Norovirus test

## 13. Quality control

Only two bands maximum must appear. **If the control band (C) is missing, the test is invalid and cannot be evaluated!** The blue control band confirms that the sample and the reagents have been added correctly, the reagents were active when the test was carried out and the sample has migrated properly through the test

membrane. A colourless background in the result zone (reaction window) acts as an internal negative control. If the test has been carried out correctly and the reagents have functioned properly, the background will be colourless in order to provide a detectable result.

# RIDA®QUICK Norovirus (N 1403)

## 14. Limitations of the method

The RIDA®QUICK Norovirus Test is used to determine Genogroup I and Genogroup II Noroviruses and confirms their presence in stool samples providing the virus loading does not fall below the detection limit of the test system. It is crucial that the stool samples are taken from the patient suffering from gastroenteritis during the phase when they are showing the symptoms. The intensity of the blue colouration of the test band does not reflect the severity of the clinical symptoms but is merely an indication of the quantity of virus antigens in the sample. The blue colouration can vary from weak to very strong.

A positive result does not rule out the presence of other pathogens. The presence of two or more pathogens simultaneously in the sample is also possible and can

only be clarified by differential diagnosis. The symptoms of double and multi-infections may be more severe than those of single causes.

A negative result does not rule out the presence of Noroviruses and their being the cause of the gastroenteritis. This may happen for the following reasons: intermittent virus elimination, the virus loading in the sample is too low, the sample was collected at the wrong time or the transport or storage conditions were not appropriate for the sample. A new sample must therefore be tested when infection with Noroviruses is suspected.

## 15. Expected values

Noroviruses cause 70 % of all non-bacterial gastroenteritis epidemics worldwide and are therefore the most frequent and are responsible for approximately 50 % enteritis outbreaks. As well as their frequent occurrence in communal facilities, they are often the cause of sporadic diarrhoeas. After Rotaviruses, they are the most frequent cause of gastroenteritis in children and affect the adults as well. Since Noroviruses vary genetically, once infected, people are often not permanently immune to re-infection.

Thus, positive test results vary in relation to the prevalence of Noroviruses within a certain population, the genotype circulating there and the environment of the patients. The way the sample was collected, transported and stored also affects the positive test result rate.

## 16. Performance characteristics

In a validation study, RIDA®QUICK Norovirus Test was tested with 113 stool samples. The samples came from patients suffering from gastroenteritis, both from sporadic and outbreak events, and were collected during the winter of 2007/2008 in the metropolitan area of Hamburg, freshly analysed using RT-PCR in the hygiene institute there and then frozen. In July 2008,

the samples were thawed and tested with the RIDA®QUICK Norovirus Test at the same institute along with the commercially available RIDASCREEN® Norovirus Elisa for comparison. The results are shown in **Table 1**.

Table 1: Correlation of the RIDA®QUICK Norovirus quick test with the RT-PCR and RIDASCREEN® Norovirus Elisa

		RIDA®QUICK		RIDASCREEN®	
		Panel WS 2007/8		Panel WS 2007/8	
		+	-	+	-
RT - PCR	+	49	8	43	13
	-	3	53	1	55

<b>Sensitivity:</b>	86.0 %	76.8 %
<b>Specificity:</b>	94.6 %	98.2 %
<b>PPV:</b>	94.2 %	97.7 %
<b>NPV:</b>	86.9 %	80.9 %
<b>Accuracy:</b>	90.3 %	87.5 %

### 16.1 Reproducibility

Three independent laboratories tested 5 stool samples each in a quadruple determination at 2 different times on one day (intra-assay variability) and on three different days (inter-assay variability). Two negative, 2 weakly positive and one medium positive sample were used.

### 16.2 Cross reactivity

Different pathogenic organisms of the intestinal tract were tested using the RIDA®QUICK Norovirus test and showed no cross reactivity. The tests were carried out on undiluted bacterial suspensions with organism

The RIDA®QUICK Norovirus Test showed 100 % reproducibility even for the control line.

concentrations ranging from 106 to 109 organisms per ml. The results are listed in **Table 2**.

Table 2: Testing for cross reactions with different micro-organisms

Test germ	Origin	Name	Result
Citrobacter freundii		DSM 30047	negative
Enterobacter cloacae		DSM 30054	negative
Enterococcus faecalis		DSM 2570	negative
Enterococcus faecium		DSM 20477	negative
E. coli	Isolate		negative
E. coli	Isolate		negative
E. hermannii		DSM 4560	negative
Lactococcus lactis		DSM 20481	negative
Listeria innocua		DSM 20649	negative
Proteus mirabilis		DSM 788	negative
Proteus mirabilis		DSM 4479	negative
Proteus vulgaris		DSM 30119	negative
Providencia stuartii			negative
Pseudomonas aeruginosa		DSM 939	negative
Salmonella Agona	Isolate		negative
Salmonella Choleraesuis		DSM 4224	negative
Salmonella Infantis	Isolate		negative
Salmonella Ohio	Isolate		negative
Salmonella Typhimurium	Isolate	DSM 554	negative
Staphylococcus aureus		DSM 20372	negative
Streptococcus agalactiae	Isolate	DSM 2134	negative

Test germ	Origin	Name	Result
Streptococcus dysgalactiae	Isolate	DSM 20662	negative
E. coli (O157:H-)	Isolate		negative
E. coli (O116:H21)	Isolate		negative
E. coli (O111:H-)	Isolate		negative
E. coli (O22:H8)	Isolate		negative
E. coli (O26:H11)	Isolate		negative
Candida albicans		ATCC 10231	negative
Salmonella enteritidis		DSM 9898	negative
Morganella morganii		DSM 6675	negative
Shigella boydii		DSMZ 7532	negative
Listeria grayi		DSMZ 20596	negative
Listeria denitrificans		DSMZ 20603	negative
Listeria grayi		DSMZ 20601	negative
Listeria seelgeri		DSMZ 20751	negative
Listeria welshimera		DSMZ 20650	negative
Listeria ivanovii subsp. Ivanovii		DSMZ 20750	negative
Listeria monocytogenes SG 1/2a		DSMZ 20600	negative
Listeria monocytogenes SG 1		ATCC 19111	negative
Listeria monocytogenes SG 4b		ATCC 19115	negative
Enterobacter sakasaki			negative

### 16.3 Interfering substances

The following substances had no effect on the test results when they were mixed with the stool samples at the specified concentrations:

- Barium sulphate (5 % w/w),
- Loperamide (immodium 5 % w/w),
- Peptobismol (5 % v/w),
- Mucin (5 % w/w),

- Sweetener (5 % v/w),
- Human blood (5 % v/w),
- Stearic acid/palmitic acid (mixed 1:1, 40 % w/w),
- Metronidazole (0,5) injection solution (5 % v/w),
- Diclofenac (0,00263 % v/w).

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**R-Biopharm AG**